

PCcaps™

Capsules for Pre-Clinical Studies

Reaching the Clinical Trial stage is never an easy process.

You need all the help you can get.

PCcaps capsules provide a simple way for Pre-Clinical oral administration of pharmaceutical actives and formulations.

With PCcaps, you can:

- Avoid the use of solvents that may be incompatible with the active, or that cause toxicity
- Prevent regurgitation and upper GI tract irritation
- Deliver exact dose to rodents without using any excipients
- Administer early phase formulations of solids (powder/granules), semi-solids (paste), or non-aqueous liquids
- Perform controlled-release studies by enteric coating of PCcaps



NEW AND IMPROVED FEATURES

- Capsules have locking rings, for improved closure
- Kit now includes a stainless steel tamper for filling powder with poor flow properties
- Dosing syringe is more flexible to provide smooth insertion into GI tract
- A CD instruction video is provided

Small size is ideal for rats and guinea pigs. Not suitable for mice or rodents below 150g weight.



PCcaps SPECIFICATIONS

Capacity	21.8 µl
Average external diameter	Body: 2.35 mm Cap: 2.64 mm
Overall Closed Length	7.2 mm
Average weight	0.01 g
Color	White Opaque

CAPSUGEL®
Quality
People and Products Working Together™

PCcaps™ KIT



- Two boxes of 50 PCcaps capsules, dosing syringe, funnel, stand and tamper




PCcaps™ REFILL

- One box of 50 PCcaps capsules



New PCcaps™ MULTIPLE FILLER:

- more convenient way to prepare multiple doses



CAPSUGEL®
Quality
People and Products Working Together™

www.capsugel.com

Capsugel Library

BAS 180 E 2003

Americas

Greenwood, USA
Capsugel Americas
535 North Emerald Road
Greenwood,
South Carolina 29646
USA
Tel: +1-864-223-2270
+1-800-845-6973
Fax: +1-864-942-3848

Puebla, Mexico

Capsugel
Av. Division Del Norte 3443
Col. San Pablo Tepetlapa
04620 Mexico D.F.
Tel: 001-800-010-0463
+52 555 326 8470
Fax: +52 555 652 2122

Rio de Janeiro, Brazil

Capsugel
Laboratórios Pfizer Ltda.
Av. das Américas,
nº 8.445 / 1.202-1.207
Barra da Tijuca
Rio de Janeiro - RJ
Cep : 22793-081 - Brazil
Tel: +55-21- 2432-5264 /
+0800704-2445
Fax: +55-21-2432-5270

Australia

Sydney, Australia
Pfizer Pty. Ltd
Capsugel Division
32, Cawarra Road
Caringbah, Sydney
NSW 2229 Australia
Tel: +61-2-9710-6642
Fax: +61-2-9710-6283

Asia

Bangkok, Thailand
Capsule (Thailand) Co., Ltd.
1/82 M005, Rojana
Industrial Park
Rojana Road,
Tambol Khanham
Amphur U-Thai,
Ayutthaya 13210
Thailand
Tel: +6635-334000
Fax: +6635-227088

Jakarta, Indonesia

Capsugel
Jl. Raya Bogor Km. 42
P.O. Box 15/CBI
Cibinong 16916,
Jawa Barta
Indonesia
Tel: +62-21-875-2226
Fax: +62-21-875-2227

Sagamihara, Japan

Capsugel
3-36 Minamihashimoto
4-Chome Sagamihara City
Kanagawa Prefecture
Japan 229-1133
Tel: +81-42-700-6700
Fax: +81-42-700-6719

Suzhou, China

Suzhou Capsugel Ltd.
JinJihu Road
Suzhou Industrial Park
Suzhou 215021
Jiangsu, PR China
Tel: +86-512-6761-6666
Fax: +86-512-6761-8807

Europe

Bornem, Belgium
Capsugel
Rijksweg 11
2880 Bornem
Belgium
Tel: +32-3-890-05-11
Fax: +32-3-889-26-22

Colmar, France

Capsugel
10, rue Timken
68027 Colmar Cedex
France
Tel: +33-3-89-20-57-09
Fax: +33-3-89-41-48-11

GMCN 200305004